

MC-H05



Overall approx dimension: 2070mm (L) x 850mm (W) X 550mm to 850mm (H). Three sectiontop. Backrest adjustable on Screw mechanism. Leg sectioncan be telescoped under the backrest for lithotomy positions. Removable Stainless steel head & foot boards with laminated panels. Foot support for high/low chair positions. Middle & Leg section tops are removable for cleaning. A pair of upholstered mild steel lithotomy Crutches mounted on SS rods. Complete with SS side railings, SS tray & SS telescopic IV pole .Height adjustable by foot operated hydraulic pump. Handle wound trendelenburg position. Three section mattress. Synthetic body castors of 125mm wheel dia, two with brake and two without brake. Pre –treated & powder coated finish. Supplied in semi knocked down condition. ONLY BED is part of the package, other items can be bought separately.



Main Features

- Three sectiontop. Backrest adjustable on Screw mechanism. Leg sectioncan be telescoped under the backrest for lithotomy positions. Removable Stainless steel head & foot boards with laminated panels. Foot support for high/low chair positions ,
- Middle & Leg section tops are removable for cleaning. A pair of upholstered mild steel lithotomy Crutches mounted on SS rods. ,
- Complete with SS side railings, SS tray & SS telescopic IV pole .Height adjustable by foot operated hydraulic pump. Handle wound trendelenburg position. ,
- Synthetic body castors of 125mm wheel dia, two with brake and two without brake. Pre –treated & powder coated finish. Supplied in semi knocked down condition.

Technical Parameter

Dimension	2070mm(L)x850mm(W)x550mm to 850mm(H)0-60°
Position	Backrest/TB/RTB/Height
Frame Work	MS Framework
Accessories	Lithotomy Crutches,IV Pole , Stainless Steel Tray, Side railing
TOP	Laminated Top
Colour	Broken White

