

**medik**

Tel: +86 15250358182

E-mail: info@medik-medical.com

Zhangjiagang  
Medi Medical Equipment  
Co., Ltd

# MK-A03

## Patient attendant chair



MK-A03 patient attendant chair in hospital Can pull out to be a bed, which saves space and more comfortable for the caregiver. Ergonomic non-slip handrail make caregiver more comfortable when seated



### SEAT EXTENSION DESIGN

Seat can be expanded into 190cm recliner, more comfortable for caregiver



### SIDE V-TYPE CLAMP

V-type clamp, self-locking system, 4pcs small wheels leave the ground after being locked



### CASTOR

Equipped with 4pcs small casters for easy extension of the chair



### ARMREST

Non-slip armrests make caregiver more comfortable when seated



### LIFTING ROPE

Lifting rope can be easily changed from unfolded state to seat state, which is safer



### LEATHER

Flame retardant, mildew resistance, bending resistance, 100,000 times wear resistance.

### Technical parameter

Model	MK-A03
Folding Size	65*62*42cm
Unfolding Size	190*62*58cm
Packing Size	118.5*61.5*22cm
Net weight	22.5kg

### Main features

- \* Suitable for high-class patient room or office
- \* Folded to be a chair, unfolded to be a bed
- \* Frame available, PVC and high quality sponge applied for chair surface
- \* Rubber casters

