

medik

Tel: +86 15250358182

E-mail: info@medik-medical.com

Zhangjiagang
Medi Medical Equipment
Co., Ltd

MK-A04

Accompany chair



MK-A04 accompany chair is for the patient's attendant in hospital, it can be fold as a chair in day time, while unfold as a bed in night, which saves space and is safe and reliable. The elegant arc-shaped armrests are ergonomically designed, comfortable and beautiful, and flexible to use.



SEAT EXTENSION

seat can be expanded into 204cm recliner, with small wheels at the bottom to easy to pull and move the seat



BUCKLE STRUCTURE

MK-A04 accompany chair can be automatically locked when it returns to the seat position to prevent it from shaking



CASTOR

Equipped with 4pcs small casters for easy extension of the chair



INSIDE EXTENSION SPRING

spring structure makes it easier and more convenient for accompany chair to change back to the seat position



WOODEN ARMREST

wooden armrest offer deluxe and gentle touch sense

Technical parameter

Folding Size	80*71*43cm
Unfolding Size	204*71*43cm
Backrest height	98cm
Packing Size	77*60*34cm/74*74*17cm
Net weight	36kg

Main features

- * Accompany chair can be open as a simple hospital bed, it can be folded into chair.
- * The structure is made of powder coating steel
- * Covered by High quality PVC artificial leather, inside is 50mm high density foam. Color is optional
- * The accompany chair is equipped with Inside extension spring
- * wooden armrest offer deluxe and gentle touch sense
- * Two pieces 3 inch noiseless wheels, four small wheels

