

# MK-A08

Accompany chair



*Better care  
Better comfort!*

# MK-A08

# ACCOMPANY CHAIR

## DESCRIPTION

Modern fashion design, excellent finishing details such as wooden hand rail, the MK-A08 acts as a comfortable and luxurious sofa chair in the daytime, and expandable as a sleep chair during night when it is needed.

## FEATURES

- \* Accompany chair can be open as a simple hospital bed, it can be folded into chair.
- \* The structure is made of powder coating steel
- \* Covered by High quality PVC artificial leather, inside is 50mm high density foam. Color is optional
- \* The accompany chair is equipped with Inside extension spring wooden armrest offer deluxe and gentle touch sense
- \* Two pieces 3 inch noiseless wheels, four small wheels

## STANDARD SPECIFICATION

Folding Size	82*86*42cm
Unfolding Size	198*86*42cm
Backrest height	89cm



**Seat Extension**

chair can be expanded into 198cm recliner, with small wheels at the bottom to easy to pull and move the chair.



**Handle**

A metal handle in front of the seat allows the chair to be extended into a recliner.



**Armrest**

Curved wooden armrests are more ergonomically designed, more comfortable.



**Buckle Structure**

There are 2 buckle rails on the side of the chair, which are used to fix the chair when folding it and prevent it from sliding.

## ABOUT US

Zhangjiagang Medi Medical Equipment Co.,Ltd, founded in 1995, with more than 20 years experience in developing hospital furniture, is one of the largest manufacturer of hospital bed, patient stretcher, medical cart and other related hospital furniture in China or even Asia.

Through scientific management complying with ISO9001&ISO13485, by adopting advanced facility like Germany TRUMPF CNC, Japan OTC Robot Welding, Auto-power coating line with Swiss GEMA gun , Medi focuses on providing brand quality product with CE approval.

Headquartered in Zhangjiagang nearby Shanghai, with subsidiary oversea offices in Hungary, Dubai, Thailand, New Zealand, Kenya, Colombia, Medi acting professionally and swiftly for both pre-sales and after-sales.The products are widely accepted by global clients, such as UK RIE hospital, UNIVERZITETNI KLINICNI CENTER, Kokilaben Dhirubhai Ambani Hospital, Singapore National University Hospital etc.



YOU TUBE



MEDIK

**Pls scan the QR code for more info**

**info@medik-medical.com | +86 139 6223 1292**

Medik reserves the right to make changes without notice in design, specifications and models. The only warranty Medik makes is the express written warranty extended on the sale or rental of its products.

© 2022 Medik Services, Inc. ALL RIGHTS RESERVED. APR53703 rev 1 01-JAN-2022 ENG – US