

# MK-B03

Hospital Crib for Newborn



*Better care  
Better comfort!*

# MK-B03

# HOSPITAL CRIB FOR NEWBORN

## DESCRIPTION

MK-B03 is a sturdy and corrosion-resistant baby crib with complete steel body and environmental protection coating. ABS bassinet has a high transparent durable body frame with great transparent visibility. With angle adjusting mechanism for needs of medical care, simple, durable and practical as a baby crib.

## FEATURES


- Steel frame with environmental protection coating ,antibacterial, corrosion resistant and durable.
- Removable ABS bassinet,good transparency, durable and safe.
- Trendelenburg position adjusted by gas spring.
- Removable utility basket for easy storage.
- Four silent castors ,dia 100mm ,with individual brakes.
- Waterproof mattress ,breathable, antibacterial and non-toxic.

## ACCESSORIES

- 1 pc of bassinet
- 1 set of waterproof mattress
- Angle adjustment handle 1 set
- 4 sets of brake casters
- 1 set of sundries net basket
- Mosquito net (optional)

## STANDARD SPECIFICATION

Overall size	L 830mm*W 500mm*H 980mm
bassinet size	L 830mm*W 500mm*H 255mm
Bassinet Safe loading	15kg
Basket Safe loading	5kg
Castor	Dia 100mm ,with individual brakes



**Angle Adjustable**

Has trendelenburg position which is adjusted by gas spring ,and the crib is removable.



**Utility Basket**

Removable utility basket on the bottom for easy storage,save space.



**Castor**

Four silent castors with individual brakes,easy to move and stand still.



**Mattress**

Thick waterproof mattress with removable cover, breathable, antibacterial and non-toxic.

## ABOUT US

Zhangjiagang Medi Medical Equipment Co.,Ltd, founded in 1995, with more than 20 years experience in developing hospital furniture, is one of the largest manufacturer of hospital bed, patient stretcher, medical cart and other related hospital furniture in China or even Asia.

Through scientific management complying with ISO9001&ISO13485, by adopting advanced facility like Germany TRUMPF CNC, Japan OTC Robot Welding, Auto-power coating line with Swiss GEMA gun , Medi focuses on providing brand quality product with CE approval.

Headquartered in Zhangjiagang nearby Shanghai, with subsidiary oversea offices in Hungary, Dubai, Thailand, New Zealand, Kenya, Colombia, Medi acting professionally and swiftly for both pre-sales and after-sales.The products are widely accepted by global clients, such as UK RIE hospital, UNIVERZITETNI KLINICNI CENTER, Kokilaben Dhirubhai Ambani Hospital, Singapore National University Hospital etc.



YOU TUBE



MEDIK

**Pls scan the QR code for more info**

**info@medik-medical.com | +86 139 6223 1292**

Medik reserves the right to make changes without notice in design, specifications and models. The only warranty Medik makes is the express written warranty extended on the sale or rental of its products.

© 2022 Medik Services, Inc. ALL RIGHTS RESERVED. APR53703 rev 1 01-JAN-2022 ENG – US