



**MODEL : MK-B07**

BABY CRIB

# BABY CRIB (MK-B07)

## MAIN FEATURES

- Removable molded Plastic Crib (without joints).
- With Trendelenburg (TR) Position. (Inclination up to 15°)
- Framework made of stainless steel.
- 1 S.S. shelf at the bottom.
- 4 swivel castor (2 with brakes)
- With a 3cm mattress, easy to clean.

## TECHNICAL PARAMETER

ITEM	PARAMETER
Overall size	820*500*860mm
Trendelenburg	0~15°
Caster	Dia75mm
Loading capacity	20kg

MK-B07 is a sturdy and corrosion-resistant baby crib with complete stainless body. Acrylic bassinet has a high transparent body with durable frame. With angle adjusting mechanism for needs of medical care, simple, durable and practical as a baby crib.




## REVERSE TRENDELENBURG ADJUST

Three angle adjusting options via simple and sturdy bolt mechanism to prevent the baby from spitting up.



### Sleeping basin and mattress

Basin use high quality ABS material and mattress is waterproof and breathable.



### BASKET

Basket under the bed frame, can store some items



### CASTOR

Dia100mm high-quality castor and two with brakes, castors are stability and safety



## Zhangjiagang Medi Medical Equipment Co., Ltd.

---

Tel : +86 13962231292 / Fax : +86 51258688450

E-mail : info@medik-medical.com

Address : No.3 Jinxiu Road,Jinfeng Industrial Zone,Jinfeng Town,Zhangjiagang,Ji-  
angsu Province,China

---



You Tube



Medik