

medik
WWW.MEDIK-MEDICAL.COM

EMERGENCY CART

MODEL:MK-C01



CONVENIENT | SECURITY | STRONG | COMFORTABLE

PRODUCT DETAILS

BETTER CARE BETTER COMFORT



» Top board

With stainless steel guard railing for better protection

» Extended Work Surface

Easy for care giver writing



» Drawer

In total 5 drawers
Inner with partitions can be well organized

» Breakaway Lock

Disposable sealing central locking system



» Stainless Steel Basket

» Sharp Container

PRODUCT DETAILS

BETTER CARE BETTER COMFORT



» IV Pole

Height adjustable IV Pole with 5 hooks

» Defibrillator Shelf

Defibrillator Shelf with multi position adjustment



» Waste Bin

Two waste bins with lid



» Power Socket

» Oxygen Tank Holder



» 4" Caster

Noiseless casters with heavy duty bearing design, 2 with brake



» TECHNICAL PARAMETER

ITEM	PARAMETER
Model	MK-C01
Overall size (L*W*H)	660*480*1100mm
Material	Steel / Aluminum / Polymer
Color	White / Red / Blue
Castor	Dia100mm
G.W	60KG
Package size	790*530*1025mm

» PRODUCT SHOW

MAIN FEATURES

- Structure: Four anti-scratch aluminum alloy columns with ABS body
- ABS plastic top board with stainless steel side rails around for protection
- With disposable sharp container, dust basket, sliding writing shelf, utility container
- Four 100mm castors and two with brakes
- Disposable sealing central locking system

