

MK-C02B

Anesthesia Cart



*Better care
Better comfort!*

MK-C02B ANESTHESIA CART

DESCRIPTION

MK-C02B anesthesia cart are with features of stability, smooth mobility, easy transportation within the hospital. They provide ample and effective storage space, with options for specialized accessories that help organize essential items such as sedatives, narcotics, and anesthetics.

FEATURES

- Epoxy Coated Steel Construction
- Convex top board with stainless steel guard railing for better protection
- 4" swivel caster wheels,two with brake
- Protective ABS bumpers on each corner
- 6 drawers
- 10 compartment tilt bin organizer

STANDARD SPECIFICATION

Model	MK-C02B
Cart Overall Size	708*560*1270mm
Tilting box	115*98*130-170mm
Structure	Aluminum alloy frame and cold steel sprayed with plastic
Castor	Dia 100mm



Drawer divide

Each drawer can be divided by different dividers, according to your needs



Tilting box

The tilting box can effectively classify and place different anesthesia supplies



Sharps Box

Sharps Box are safer for disposing of medical syringes, needles and blades.



Waste bin

Standard with flip-top waste bin, optional foot pedal waste bin.



Glove Storage Box

The glove storage box is convenient and practical, improving the work efficiency of medical staff



Mobile holder

Convenient for medical staff to adjust positions and configure relevant containers according to needs



Castor

Luxury silent 4 inch casters, 2 with brake



Central lock

One key controls the locking of multiple drawers, making it convenient for professionals to manage and use them.

ABOUT US

Zhangjiagang Medi Medical Equipment Co.,Ltd, founded in 1995, with more than 20 years experience in developing hospital furniture, is one of the largest manufacturer of hospital bed, patient stretcher, medical cart and other related hospital furniture in China or even Asia.

Through scientific management complying with ISO9001&ISO13485, by adopting advanced facility like Germany TRUMPF CNC, Japan OTC Robot Welding, Auto-power coating line with Swiss GEMA gun , Medi focuses on providing brand quality product with CE approval.

Headquartered in Zhangjiagang nearby Shanghai, with subsidiary oversea offices in Hungary, Dubai, Thailand, New Zealand, Kenya, Colombia, Medi acting professionally and swiftly for both pre-sales and after-sales.The products are widely accepted by global clients, such as UK RIE hospital, UNIVERZITETNI KLINICNI CENTER, Kokilaben Dhirubhai Ambani Hospital, Singapore National University Hospital etc.



YOU TUBE



MEDIK

Pls scan the QR code for more info

info@medik-medical.com | +86 139 6223 1292

Medik reserves the right to make changes without notice in design, specifications and models. The only warranty Medik makes is the express written warranty extended on the sale or rental of its products.

© 2022 Medik Services, Inc. ALL RIGHTS RESERVED. APR53703 rev 1 01-JAN-2022 ENG – US