

MK-C03A

Medication Cart



*Better care
Better comfort!*

MK-C03A

MEDICATION CART

DESCRIPTION

Our MK-C03A Medication Cart focuses on practical and functional design to meet the needs of different departments. The flexible partition system, removable drawer trays and double-door storage cabinets provide ample storage space for a variety of medical equipment and supplies.

FEATURES

- Made of aluminum alloy frame and cold steel sprayed with plastic.
- The top can be opened and there's a storage space inner with dividers, which can be used to store large items.
- One ABS push handle, two garbage cans with different colors, and one sharp box on the left.
- Four 100mm diameter luxury castors and two of them with brakes.
- Front equipped with two layers of drawers and center lock system and double doors cabinet, all with inner dividers.

TECHNICAL PARAMETER

Model	MK-C03A
Overall size (L*W*H)	620*420*900mm
Castor	Dia100mm
Optional	Sharps box basket, removable IV pole, CPR board
Packing size	790*540*840mm
G.W	31.2kg



ABS handle

ABS handle makes it easier for medical staff to move the cart.



Hidden storage space

The top can be opened and there is a storage space with dividers inside for storing large items.



Waste bin

Standard with flip-top waste bin, optional foot pedal waste bin.



Castors

Luxury silent 4 inch casters, 2 with brake.



Column control

The main controller is installed on the column and contains removable battery, height controller and emergency stop button.



Spreader Bar

4 hooks spreader Bar can firm and stable lifting sling.



ABOUT US

Zhangjiagang Medi Medical Equipment Co.,Ltd, founded in 1995, with more than 20 years experience in developing hospital furniture, is one of the largest manufacturer of hospital bed, patient stretcher, medical cart and other related hospital furniture in China or even Asia.

Through scientific management complying with ISO9001&ISO13485, by adopting advanced facility like Germany TRUMPF CNC, Japan OTC Robot Welding, Auto-power coating line with Swiss GEMA gun , Medi focuses on providing brand quality product with CE approval.

Headquartered in Zhangjiagang nearby Shanghai, with subsidiary oversea offices in Hungary, Dubai, Thailand, New Zealand, Kenya, Colombia, Medi acting professionally and swiftly for both pre-sales and after-sales.The products are widely accepted by global clients, such as UK RIE hospital, UNIVERZITETNI KLINICNI CENTER, Kokilaben Dhirubhai Ambani Hospital, Singapore National University Hospital etc.



YOU TUBE



MEDIK

Pls scan the QR code for more info

info@medik-medical.com | +86 139 6223 1292

Medik reserves the right to make changes without notice in design, specifications and models. The only warranty Medik makes is the express written warranty extended on the sale or rental of its products.

© 2022 Medik Services, Inc. ALL RIGHTS RESERVED. APR53703 rev 1 01-JAN-2022 ENG – US