

# MK-C03B

Medication Cart



*Better care  
Better comfort!*

# MK-C03B

## MEDICATION CART

### DESCRIPTION

MK-C03B Medication Cart focuses on practical and functional design to meet the needs of different departments. Specially the top sliding door make the workflow very convenient and easy to access the urgent medical materials. The flexible partition system, double-door storage cabinets provide ample storage space for a variety of medical equipment and supplies.

### FEATURES

- Made of aluminum alloy frame and cold steel sprayed with plastic.
- The top can be opened and there's a storage space inner with dividers, which can be used to store large items.
- Four 100mm diameter luxury castors and two of them with brakes.
- Front equipped with two layers of drawers and center lock system and double doors cabinet, all with inner dividers.

### TECHNICAL PARAMETE

Model	MK-C03B
Overall size (L*W*H)	620*420*900mm
Castor	Dia100mm
Optional	Sharps box basket, removable IV pole
Packing size	790*540*840mm
G.W	31.2kg



### Storage Space And Push Handle

Detachable storage basket and ABS push handle on the right side of the cart .



### Hidden Storage Space

The top can be opened and there is a storage space with dividers inside for storing large items.



### Waste Bin

Standard with flip-top waste bin, optional foot pedal waste bin.



### Castor

Luxury silent 4 inch casters, 2 with brake.



### Open Top

The top plate is equipped with aluminum alloy guardrails on three sides to prevent items from falling. The key lock controls the state of the sliding top plate, which is convenient for medical staff to effectively manage medical items.



### Double Doors Cabinet

A partition inside the door can be used to store small instruments.

## ABOUT US

Zhangjiagang Medi Medical Equipment Co.,Ltd, founded in 1995, with more than 20 years experience in developing hospital furniture, is one of the largest manufacturer of hospital bed, patient stretcher, medical cart and other related hospital furniture in China or even Asia.

Through scientific management complying with ISO9001&ISO13485, by adopting advanced facility like Germany TRUMPF CNC, Japan OTC Robot Welding, Auto-power coating line with Swiss GEMA gun , Medi focuses on providing brand quality product with CE approval.

Headquartered in Zhangjiagang nearby Shanghai, with subsidiary oversea offices in Hungary, Dubai, Thailand, New Zealand, Kenya, Colombia, Medi acting professionally and swiftly for both pre-sales and after-sales.The products are widely accepted by global clients, such as UK RIE hospital, UNIVERZITETNI KLINICNI CENTER, Kokilaben Dhirubhai Ambani Hospital, Singapore National University Hospital etc.



YOU TUBE



MEDIK

**Pls scan the QR code for more info**

**info@medik-medical.com | +86 139 6223 1292**

Medik reserves the right to make changes without notice in design, specifications and models. The only warranty Medik makes is the express written warranty extended on the sale or rental of its products.

© 2022 Medik Services, Inc. ALL RIGHTS RESERVED. APR53703 rev 1 01-JAN-2022 ENG – US