

MK-C08

Medical Trolley



*Better care
Better comfort!*

MK-C08

MEDICA TROLLY

DESCRIPTION

MK-C08 not only has the aesthetic appearance, but also has a more practical design on function. It will bring the efficient working experience to employees with stable and smooth use, convenient operation and perfect appearance.

FEATURES

- Made of aluminum alloy frame and cold steel sprayed with plastic.
- Front equipped with four layers of drawers and has center lock system,3 middle and 1big ,all have dividers inner.
- Top can be opened and there's a storage space inner with diviers,which can be used to store large items.
- One ABS push handle and two detachable storage boxes on the right of the cart ,two garbage cans with different colors and one sharp box on the left of the cart.
- Four 100mm diameter luxury castors and two of them with brake.

STANDARD SPECIFICATION

Overall Size	660*480*960MM
Structure	aluminum alloy frame and cold steel sprayed with plastic
Castor	Dia 100mm



Open top

Three sides of top plate has aluminum handles, can be both easy to push the cart and prevent stuffs from dropping down. Sliding open top plate can let medical staff to access the medicine stored underneath very easily and immediately. The storage space on the top with dividers inner which can be separated freely while it can be used for the storage of ambulance bag and Infusion bottle either if without dividers.



Locking system

Front key can lock the top mesa, the center key lock on this side is used for front drawers.




Drawer

Three middle drawers (size: 470*365*110mm) and one big drawer (size: 470*365*145mm) on the front, both of them have an aluminum handle, there are dividers and sealing slot-type signage in the drawers which can be separated freely and easily to recognize.



Garbage can

Two garbage cans with different colors on the left side of the cart which can separate the garbage, also there is a sharp box there to collect syringes.



Castor

There are four 100mm diameter luxury castors on the bottom and two of them with a brake.



Storage space and push handle

Two detachable storage spaces and ABS push handle on the right side of the cart.



Anti-collision

Anti-collision design around the cart.

ABOUT US

Zhangjiagang Medi Medical Equipment Co.,Ltd, founded in 1995, with more than 20 years experience in developing hospital furniture, is one of the largest manufacturer of hospital bed, patient stretcher, medical cart and other related hospital furniture in China or even Asia.

Through scientific management complying with ISO9001&ISO13485, by adopting advanced facility like Germany TRUMPF CNC, Japan OTC Robot Welding, Auto-power coating line with Swiss GEMA gun , Medi focuses on providing brand quality product with CE approval.

Headquartered in Zhangjiagang nearby Shanghai, with subsidiary oversea offices in Hungary, Dubai, Thailand, New Zealand, Kenya, Colombia, Medi acting professionally and swiftly for both pre-sales and after-sales.The products are widely accepted by global clients, such as UK RIE hospital, UNIVERZITETNI KLINICNI CENTER, Kokilaben Dhirubhai Ambani Hospital, Singapore National University Hospital etc.



YOU TUBE



MEDIK

Pls scan the QR code for more info

info@medik-medical.com | +86 139 6223 1292

Medik reserves the right to make changes without notice in design, specifications and models. The only warranty Medik makes is the express written warranty extended on the sale or rental of its products.

© 2022 Medik Services, Inc. ALL RIGHTS RESERVED. APR53703 rev 1 01-JAN-2022 ENG – US