

# MK-P01

Emergency carts



*Better care  
Better comfort!*

# MK-P01

# Emergency carts

## DESCRIPTION

With ergonomic design, fascinating appearance, the MK-P01 emergency cart is with very identical features. The cart is manufactured with stable structure for harsh manipulation, all excellent features with optional accessories provide the medical staff a lovely working platform.

## FEATURES

- Division for the drawer
- Holding tray for equipment as defibrillators
- Oxygen tank holder
- CPR board
- Serum holder
- Slide-out work surface
- Precision-welded steel frame with durable, high-impact panels
- Delivered fully assembled, no tools required

## STANDARD SPECIFICATION

### CART SPECIFICATION

	Height	Depth	Width
Large	110cm/43"	61cm/24"	79cm/31"
Standard	100cm/39.5"	61cm/24"	79cm/31"
Small	91cm/36"	61cm/24"	79cm/31"

### DRAWER SPECIFICATION

	Height	Depth	Width
Large	7cm/3"	45cm/17.5"	79cm/31"
Standard	16cm/6"	45cm/17.5"	79cm/31"
Small	25cm/10"	45cm/17.5"	79cm/31"



**Drawer divider**

Each drawer can be divided by different dividers, according to your needs.



**Defibrillator shelf**

The shelf can hold emergency equipment such as a defibrillator.



**Extended work surface**

The pull-out extension surface is convenient for medical staff to write on it.



**Oxygen cylinder holder**

Oxygen cylinders can be placed for emergency use.



**CPR board**

Size 475\*595\*8.5mm



**Waste bin**

Standard with flip-top waste bin, optional foot pedal waste bin.



**IV pole**

2 hooks, height adjustable.



**Castor**

Luxury silent 5 inch casters, 2 with brake.



**Suction shelf(option)**

## ABOUT US

Zhangjiagang Medi Medical Equipment Co.,ltd, founded in 1995, with more than 20 years experience in developing hospital furniture, is one of the largest manufacturer of hospital bed, patient stretcher, medical cart and other related hospital furniture in China or even Asia.

Through scientific management complying with ISO9001&ISO13485, by adopting advanced facility like Germany TRUMPF CNC, Japan OTC Robot Welding, Auto-power coating line with Swiss GEMA gun , Medi focuses on providing brand quality product with CE approval.

Headquartered in Zhangjiagang nearby Shanghai, with subsidiary oversea offices in Hungary, Dubai, Thailand, New Zealand, Kenya, Colombia, Medi acting professionally and swiftly for both pre-sales and after-sales.The products are widely accepted by global clients, such as UK RIE hospital, UNIVERZITETNI KLINICNI CENTER, Kokilaben Dhirubhai Ambani Hospital, Singapore National University Hospital etc.



YOU TUBE



MEDIK

**Pls scan the QR code for more info**

**info@medik-medical.com | +86 139 6223 1292**

Medik reserves the right to make changes without notice in design, specifications and models. The only warranty Medik makes is the express written warranty extended on the sale or rental of its products.

© 2022 Medik Services, Inc. ALL RIGHTS RESERVED. APR53703 rev 1 01-JAN-2022 ENG – US