

MK-P07

Anesthesia carts



*Better care
Better comfort!*

MK-P07

Anesthesia carts

DESCRIPTION

MK-P07 is designed with a synergy of function and technology to meet the unique needs of your anesthesia department. A flexible divider system, removable drawer trays, and a host of organizer accessories neatly arrange the medications and all related supply items.

FEATURES

- Division for the drawer
- Slide-out work surface
- Multi upper tilting bin container
- Precision-welded steel frame with durable, high-impact panels
- Delivered fully assembled, no tools required

STANDARD SPECIFICATION

CART SPECIFICATION

	Height	Depth	Width
Large	110cm/43"	61cm/24"	79cm/31"
Standard	100cm/39.5"	61cm/24"	79cm/31"
Small	91cm/36"	61cm/24"	79cm/31"

DRAWER SPECIFICATION

	Height	Depth	Width
Large	25cm/10"	45cm/17.5"	55cm/21.5"
Middle	16cm/6"	45cm/17.5"	55cm/21.5"
Small	7cm/3"	45cm/17.5"	55cm/21.5"



DRAWER DIVIDER

Each drawer can be divided by different dividers, according to your needs.



TILTING BIN

The transparent bins are above the cart, and the medicines placed can be found quickly.



WASTE BIN

Standard with flip-top waste bin, optional foot pedal waste bin.



EXTENDED WORK SURFACE

The pull-out extension surface is convenient for medical staff to write on it.



CASTOR

Luxury silent 5 inch casters, 2 with brake.



DISPOSABLE SEAL(STANDARD)

ABOUT US

Zhangjiagang Medi Medical Equipment Co.,Ltd, founded in 1995, with more than 20 years experience in developing hospital furniture, is one of the largest manufacturer of hospital bed, patient stretcher, medical cart and other related hospital furniture in China or even Asia.

Through scientific management complying with ISO9001&ISO13485, by adopting advanced facility like Germany TRUMPF CNC, Japan OTC Robot Welding, Auto-power coating line with Swiss GEMA gun , Medi focuses on providing brand quality product with CE approval.

Headquartered in Zhangjiagang nearby Shanghai, with subsidiary oversea offices in Hungary, Dubai, Thailand, New Zealand, Kenya, Colombia, Medi acting professionally and swiftly for both pre-sales and after-sales. The products are widely accepted by global clients, such as UK RIE hospital, UNIVERZITETNI KLINICNI CENTER, Kokilaben Dhirubhai Ambani Hospital, Singapore National University Hospital etc.



YOU TUBE



MEDIK

Pls scan the QR code for more info

info@medik-medical.com | +86 139 6223 1292

Medik reserves the right to make changes without notice in design, specifications and models. The only warranty Medik makes is the express written warranty extended on the sale or rental of its products.

© 2022 Medik Services, Inc. ALL RIGHTS RESERVED. APR53703 rev 1 01-JAN-2022 ENG – US