

MK-P11

Medical Trolley



*Better care
Better comfort!*

MK-P11

Medical Trolley

DESCRIPTION

Not only easy manipulation, large storage, but also accessible from both sides of the cart, MK-P11 medicine cart provides the medical staff a fantastic operation experience with efficient workflow.

FEATURES

1. Made of aluminum alloy frame and ABS plastic structure, light and firm
2. Double-sided medicine storage space, big drawer on the front of the cart, the large storage space provides a great convenience for nurses.
3. Detachable storage box and two-color garbage can.
4. Center-lock system both for front and back
5. four 100mm diameter luxury castors, two of them with brake

STANDARD SPECIFICATION

| | |
|--------------|--|
| Overall Size | 750*475*950mm |
| Structure | aluminum alloy frame and ABS plastic structure |
| Castor | Dia 100mm |



Garbage can

There are two garbage can with different colors on the left of the cart which can be used for garbage sorting.



Storage space

There is a detachable storage box on the right side of the cart.



Mesa

ABS guardrail one-time molding gentle slope design, handrail rail dual-use integration design, it can prevent the goods from slipping and make the cart easy to move.



Drawer

Both front and back have five-layer drawers with smooth track have dividers and insert drawers. labels. Each side has 25 d



Locking system

Both front and back are center lock of key.



Castor

There are four 100mm diameter and luxury castors on the bottom, two of them with brake.

ABOUT US

Zhangjiagang Medi Medical Equipment Co.,Ltd, founded in 1995, with more than 20 years experience in developing hospital furniture, is one of the largest manufacturer of hospital bed, patient stretcher, medical cart and other related hospital furniture in China or even Asia.

Through scientific management complying with ISO9001&ISO13485, by adopting advanced facility like Germany TRUMPF CNC, Japan OTC Robot Welding, Auto-power coating line with Swiss GEMA gun , Medi focuses on providing brand quality product with CE approval.

Headquartered in Zhangjiagang nearby Shanghai, with subsidiary oversea offices in Hungary, Dubai, Thailand, New Zealand, Kenya, Colombia, Medi acting professionally and swiftly for both pre-sales and after-sales. The products are widely accepted by global clients, such as UK RIE hospital, UNIVERZITETNI KLINICNI CENTER, Kokilaben Dhirubhai Ambani Hospital, Singapore National University Hospital etc.



YOU TUBE



MEDIK

Pls scan the QR code for more info

info@medik-medical.com | +86 139 6223 1292

Medik reserves the right to make changes without notice in design, specifications and models. The only warranty Medik makes is the express written warranty extended on the sale or rental of its products.

© 2022 Medik Services, Inc. ALL RIGHTS RESERVED. APR53703 rev 1 01-JAN-2022 ENG – US