

MK-PC01

Mobile Medical Computer Cart



*Better care
Better comfort!*

MK-PC01

MOBILE MEDICAL COMPUTER CART

DESCRIPTION

MK-PC01 is a mobile medical cart designed for information management during dairy patient checkup, it can bring the stable and smooth sense of use to doctors with its ergonomic and humanized design and durable quality, doctors can record the data of patients anytime during the ward rounds and it can improve the work efficiency.

FEATURES

- The column is made of aluminum alloy, and the height adjustment range is 210mm.
- The car body is made of cold steel sprayed plastic, the working surface and drawers are integrally formed with ABS engineering plastics, durable and easy to clean.
- Height of the ABS side rail is 60mm, the mouse tray can be pulled out from the left and right sides, humanized design.
- Garbage can on the front of the cart ,easy to store trash.
- Working surface with locking enclosure for laptop and has a security lock.
- Four double-sided luxury silent wheel with diagonal locking which is durable and easy to move.

STANDARD SPECIFICATION

Overall size	480*520*900/938/1138mm
Height adjustment range	210mm
Keyboard drawer spec	463*200*500mm
Mouse tray spec	200*200mm
castors	Dia 100mm



Locking enclosure for laptop

The working surface has a computer locking enclosure which can fix computers of different sizes, and is equipped with a security lock.



Height adjustment function

Effortless height adjustment with foot pedal, it can stop at different height to meet different needs which can adapted to doctors of different heights and for use both in standing and sitting positions.



Working surface

The working surface and drawers are integrally formed with ABS engineering plastics which is durable and easy to clean and both has handles. The mouse tray can be pulled out from the left and right sides which can meet different usage habits of doctors.



Garbage can

Garbage can on the front of the cart which is easy to store trash.



Castors

Four double-sided luxury silent wheel with diagonal locking which can make the cart body easy to move and stand still after locked.

ABOUT US

Zhangjiagang Medi Medical Equipment Co.,Ltd, founded in 1995, with more than 20 years experience in developing hospital furniture, is one of the largest manufacturer of hospital bed, patient stretcher, medical cart and other related hospital furniture in China or even Asia.

Through scientific management complying with ISO9001&ISO13485, by adopting advanced facility like Germany TRUMPF CNC, Japan OTC Robot Welding, Auto-power coating line with Swiss GEMA gun , Medi focuses on providing brand quality product with CE approval.

Headquartered in Zhangjiagang nearby Shanghai, with subsidiary oversea offices in Hungary, Dubai, Thailand, New Zealand, Kenya, Colombia, Medi acting professionally and swiftly for both pre-sales and after-sales. The products are widely accepted by global clients, such as UK RIE hospital, UNIVERZITETNI KLINICNI CENTER, Kokilaben Dhirubhai Ambani Hospital, Singapore National University Hospital etc.



YOU TUBE



MEDIK

Pls scan the QR code for more info

info@medik-medical.com | +86 139 6223 1292

Medik reserves the right to make changes without notice in design, specifications and models. The only warranty Medik makes is the express written warranty extended on the sale or rental of its products.

© 2022 Medik Services, Inc. ALL RIGHTS RESERVED. APR53703 rev 1 01-JAN-2022 ENG – US