

MK-PC06

Multifunction Medical
Computer Cart



*Better care
Better comfort!*

MULTIFUNCTION MEDICAL COMPUTER CART

DESCRIPTION

MK-PC06 is a cart with integrated concept of medicine trolley and patient information management. With monitor stand and battery storage well built into the cart, this multi-function cart is becoming popular in hospitals and widely accepted by medical staff for improving working efficiency as well as space saving.

FEATURES

- The car body frame is made of cold steel spray plastic material, and the working surface and drawers are made of ABS engineering plastic, which is wear-resistant and easy to clean.
- Working surface has 600mm height side rail and double push handle. Drawers with central locking on both sides of the cart, the first drawer can place the keyboard and mouse, the second and the third drawer height: 80mm, inner size: 424*375*68 mm, the fourth drawer height: 120mm, inner size: 424*375*110 mm, the fifth drawer height: 200 mm, inner size: 424*375*182 mm, the sixth drawer height: 240mm, inner size: 424*375*220 mm. All drawers have dividers can be separated freely.
- Left of the cart has a basket with two 2L sharp boxes, a hand sanitizer holder and two garbage cans with different colors, size of garbage cans: 175*175*280MM.
- Right of the cart has a foldable working surface.
- The monitor stand can hang the display to rotate 360 degrees, the tilt angle can be adjusted.
- Four double-sided luxury silent wheel with diagonal locking which is durable and easy to move, anti-collision design.

SPECIFICATION

Overall size	520*510*1020mm
Monitor stand load capacity	Standard 10kg, max option 21Kg
Castors	Dia 100mm
Garbage can size	175*175*280mm
Sharp box size	2L



Monitor Stand

The cart equipped with a monitor stand to fix the monitor which can be adjusted angle and 360-degree rotation.



Working Surface

The working surface is integrally formed of ABS engineering plastics, with double handle which is convenient for pushing cart, and the ABS guardrail is 60mm height which can prevent items from slipping.



Castors

Four double-sided luxury silent wheel with diagonal locking which can make the cart body easy to move and stand still after locked, anti-collision design around which can increase safety.



Garbage Can And Storage Basket

One storage basket with two sharp boxes on the left of the cart ,next to it is a hand sanitizer holder,two garbage cans with different color which is easy to separate medical waste.



Foldable Working Surface

The right side of the cart has a foldable working surface which can increase the usable area and flexibility of the cart.



Drawers With Central Locking

Six drawers with central locking on both sides of the cart ,the first drawer is designed for the keyboard and mouse ,the rest five drawers are storage spaces with dividers which can be separated freely ,all of the drawers adopts silent tracks.Different size of drawers can meet different needs of various departments for ward rounds.

ABOUT US

Zhangjiagang Medi Medical Equipment Co.,Ltd, founded in 1995, with more than 20 years experience in developing hospital furniture, is one of the largest manufacturer of hospital bed, patient stretcher, medical cart and other related hospital furniture in China or even Asia.

Through scientific management complying with ISO9001&ISO13485, by adopting advanced facility like Germany TRUMPF CNC, Japan OTC Robot Welding, Auto-power coating line with Swiss GEMA gun , Medi focuses on providing brand quality product with CE approval.

Headquartered in Zhangjiagang nearby Shanghai, with subsidiary oversea offices in Hungary, Dubai, Thailand, New Zealand, Kenya, Colombia, Medi acting professionally and swiftly for both pre-sales and after-sales. The products are widely accepted by global clients, such as UK RIE hospital, UNIVERZITETNI KLINICNI CENTER, Kokilaben Dhirubhai Ambani Hospital, Singapore National University Hospital etc.



YOU TUBE



MEDIK

Pls scan the QR code for more info

info@medik-medical.com | +86 139 6223 1292

Medik reserves the right to make changes without notice in design, specifications and models. The only warranty Medik makes is the express written warranty extended on the sale or rental of its products.

© 2022 Medik Services, Inc. ALL RIGHTS RESERVED. APR53703 rev 1 01-JAN-2022 ENG – US