

MK-PC07

Mobile Multi-Scene
Computer Cart



*Better care
Better comfort!*

MK-PC07

MOBILE MULTI-SCENE COMPUTER CART

DESCRIPTION

MK-PC07 is a cart with concept of integrated medicine trolley and patient information management workstation. This multi-scene cart is becoming popular in hospitals and widely accepted by medical staffs for improving working efficiency as well as space saving.

FEATURES

1. Car body frame is made of SPHC steel plate with aluminum alloy, ABS engineering plastic working surface and drawers, wear-resistant and easy to clean.
2. Three small drawers & one middle & one big drawer, abundant Storage Space.
3. Storage baskets and garbage can meet daily storage needs.
4. Extension working surface tray can increase the usable area.
5. Keyboard and mouse tray can meet different usage habits.
6. The monitor stand can fix the display and rotated, angle and height can be adjusted.
7. Four double-sided luxury silent wheel with diagonal locking, durable and easy to move, anti-collision design.
8. Back drawer can store the computer lines.

STANDARD SPECIFICATION

Model	MK-PC07
Overall size	605*576*1045mm
Working surface size	L 606*W 637mm
Max load capacity	90kg
Castors	Dia 100mm
G.W	41kg



Monitor Stand

The cart equipped with a monitor stand to fix the monitor which has height and angle adjustment function and can be rotated.



Working Surface

Working surface is ABS engineering plastics, double handle is convenient for pushing and can prevent stuffs from slipping, also equipped with a soft glass, no cleaning dead angle, easy to disinfect and clean.



Drawers With Central Locking

Six drawers with central locking on the right of the cart.



Keyboard And Mouse Tray

The first layer is designed for keyboard and mouse, it can be pulled out both right and left side which can adapt to different usage habits.



Extension Tray

There is an extension tray on the right of the cart, after being pulled out, the usable area of the working surface can be increased.



Abundant Storage Space

Three small drawers, one middle and one big drawer, with dividers which can be separated freely, all of them adopt silent tracks. Different sizes of drawers can meet different needs of various departments of hospital.



Storage Basket

Two storage baskets on the right of the cart, one of them has two sharp boxes.



ABS Storage Box And PP Garbage Can

Storage box and garbage can with different colors on the left of the cart which can meet daily storage needs.



Drawer And Hook

A drawer designed to store the electric lines of the computer, and also a metal hook for hanging items.

ABOUT US

Zhangjiagang Medi Medical Equipment Co.,Ltd, founded in 1995, with more than 20 years experience in developing hospital furniture, is one of the largest manufacturer of hospital bed, patient stretcher, medical cart and other related hospital furniture in China or even Asia.

Through scientific management complying with ISO9001&ISO13485, by adopting advanced facility like Germany TRUMPF CNC, Japan OTC Robot Welding, Auto-power coating line with Swiss GEMA gun, Medi focuses on providing brand quality product with CE approval.

Headquartered in Zhangjiagang nearby Shanghai, with subsidiary oversea offices in Hungary, Dubai, Thailand, New Zealand, Kenya, Colombia, Medi acting professionally and swiftly for both pre-sales and after-sales. The products are widely accepted by global clients, such as UK RIE hospital, UNIVERZITETNI KLINICNI CENTER, Kokilaben Dhirubhai Ambani Hospital, Singapore National University Hospital etc.



YOU TUBE



MEDIK

Pls scan the QR code for more info

info@medik-medical.com | +86 139 6223 1292

Medik reserves the right to make changes without notice in design, specifications and models. The only warranty Medik makes is the express written warranty extended on the sale or rental of its products.

© 2022 Medik Services, Inc. ALL RIGHTS RESERVED. APR53703 rev 1 01-JAN-2022 ENG – US