

MK-PC08

Mobile Computer Cart



*Better care
Better comfort!*

MK-PC08

MOBILE COMPUTER CART

DESCRIPTION

MK-PC08 is a new generation of mobile medical cart designed for information management for hospital to improve the work flow efficiency, it can bring the stable and smooth sense of use to doctors with its ergonomic design and durable quality, doctors can record the data of patients anytime for dairy room checkup or in nurse station.

FEATURES

1. The car body frame is made of SPHC steel plate with aluminum alloy, the working surface and drawers are made of ABS engineering plastic , which is wear-resistant and easy to clean.
2. Keyboard and mouse tray can meet different usage habits.
- 3.The monitor stand can fix the display , can be rotated as well as the angle and height can be adjusted.
4. Four double-sided luxury silent wheel with diagonal locking which is durable and easy to move,anti-collision design.

STANDARD SPECIFICATION

Model	MK-PC08
Overall size	550*510*780~1080mm
Working surface size	L550*W510mm
Max load capacity	90kg
Castors	Dia 100mm
Keyboard tray size	L480*W130mm
G.W	41kg



Monitor Stand

The cart equipped with a monitor stand to fix the monitor which has height and angle adjustment function and can be rotated.



Working Surface

The working surface is integrally formed of ABS engineering plastics, with double handle which is convenient for pushing cart and can prevent stuffs from slipping, also equipped with a soft glass, easy to disinfect and clean.



Keyboard And Mouse Tray

The first layer is designed for keyboard and mouse, the mouse tray can be pulled out both right and left side which can adaptable to different usage habits.



Removable Storage Basket

A small storage basket is equipped together with the mouse tray which can hold small stuffs which is convenient and can be removable.



Adjustment Handle

It uses a gas spring for smooth height adjustment.



Castors

Four double-sided luxury silent wheel with diagonal locking which can make the cart body easy to move and stand still after locked, anti-collision design around which can increase safety.

ABOUT US

Zhangjiagang Medi Medical Equipment Co.,Ltd, founded in 1995, with more than 20 years experience in developing hospital furniture, is one of the largest manufacturer of hospital bed, patient stretcher, medical cart and other related hospital furniture in China or even Asia.

Through scientific management complying with ISO9001&ISO13485, by adopting advanced facility like Germany TRUMPF CNC, Japan OTC Robot Welding, Auto-power coating line with Swiss GEMA gun , Medi focuses on providing brand quality product with CE approval.

Headquartered in Zhangjiagang nearby Shanghai, with subsidiary oversea offices in Hungary, Dubai, Thailand, New Zealand, Kenya, Colombia, Medi acting professionally and swiftly for both pre-sales and after-sales.The products are widely accepted by global clients, such as UK RIE hospital, UNIVERZITETNI KLINICNI CENTER, Kokilaben Dhirubhai Ambani Hospital, Singapore National University Hospital etc.



YOU TUBE



MEDIK

Pls scan the QR code for more info

info@medik-medical.com | +86 139 6223 1292

Medik reserves the right to make changes without notice in design, specifications and models. The only warranty Medik makes is the express written warranty extended on the sale or rental of its products.

© 2022 Medik Services, Inc. ALL RIGHTS RESERVED. APR53703 rev 1 01-JAN-2022 ENG – US