

MK-PK02

pediatric emergency cart



*Better care
Better comfort!*

MK-PK02

pediatric emergency cart

DESCRIPTION

MK-PK02 is a pediatric emergency cart designed for children's treatment environment, which is more reasonable after proven clinical methodology. Color coded drawers enable doctors and clinicians to quickly identify the medicines when every seconds count.

FEATURES

- Color-coded drawers to quickly identify critical supplies when responding to emergencies.
- Slide-out work surface
- An integral divider system that can be changed and rearranged
- Keyless Access with manual relocking electric lock system
- Wide selection of accessories and organization systems
- Delivered fully assembled, no tools required

STANDARD SPECIFICATION

CART SPECIFICATION

	Height	Depth	Width
Large	110cm/43"	61cm/24"	79cm/31"
Standard	100cm/39.5"	61cm/24"	79cm/31"
Small	91cm/36"	61cm/24"	79cm/31"

DRAWER SPECIFICATION

	Height	Depth	Width
Large	7cm/3"	45cm/17.5"	79cm/31"
Standard	16cm/6"	45cm/17.5"	79cm/31"
Small	25cm/10"	45cm/17.5"	79cm/31"



Electronic lock system

The safe and convenient electronic lock can have a variety of password options, and the operation is simple and convenient. Keyless Access with manual Relocking.



Color-coded drawers

There are different color stickers on the drawer, help quickly identify critical supplies when responding to emergencies.



Dividers

Each drawer has adjustable dividers that can be adjusted as needed.



Waste container

Large capacity waste container with flap.



Lockable castors

Luxury silent 5 inch casters, 2 with brake.



Dividers

Each drawer has adjustable dividers that can be adjusted as needed.



Single Glove box

The glove carton can be stored in the box for easy placement and extraction.



Defibrillator shelf

The 90-degree rotating Defibrillator Shelf can have three gears and a strap fixing device.



Aspirator bracket

Retractable aspirator stand with straps to prevent the aspirator from falling.

ABOUT US

Zhangjiagang Medi Medical Equipment Co.,Ltd, founded in 1995, with more than 20 years experience in developing hospital furniture, is one of the largest manufacturer of hospital bed, patient stretcher, medical cart and other related hospital furniture in China or even Asia.

Through scientific management complying with ISO9001&ISO13485, by adopting advanced facility like Germany TRUMPF CNC, Japan OTC Robot Welding, Auto-power coating line with Swiss GEMA gun , Medi focuses on providing brand quality product with CE approval.

Headquartered in Zhangjiagang nearby Shanghai, with subsidiary oversea offices in Hungary, Dubai, Thailand, New Zealand, Kenya, Colombia, Medi acting professionally and swiftly for both pre-sales and after-sales. The products are widely accepted by global clients, such as UK RIE hospital, UNIVERZITETNI KLINICNI CENTER, Kokilaben Dhirubhai Ambani Hospital, Singapore National University Hospital etc.



YOU TUBE



MEDIK

Pls scan the QR code for more info

info@medik-medical.com | +86 139 6223 1292

Medik reserves the right to make changes without notice in design, specifications and models. The only warranty Medik makes is the express written warranty extended on the sale or rental of its products.

© 2022 Medik Services, Inc. ALL RIGHTS RESERVED. APR53703 rev 1 01-JAN-2022 ENG – US