



## MAIN FEATURES

Next-Level Innovation for Next-Level Care

Dual-Bar Design: Help prevent medial leg drop with this innovative safety feature

Boot Release Handle: Easily adjust leg position with a single-point latch  
Welded Pads, Removable Silicone Straps: Control potential infection with easy-to-clean welded pads and silicone straps

Patient positioning using the Lithotomy Stirrups features an extended lateral boot design that supports the patient's leg entirely, preventing superficial nerve damage and prevents unwanted leg rotation intra-operatively.

**medik**

Tel : +86 15250358182

E-mail : info@medik-medical.com



# YA-01S

## LITHOTOMY STIRRUP



## TECHNICAL PARAMETER

ITEM	PARAMETER
Cantilever Length	1000mm
Back section up/down	90° / -35°
Lift up and down	30° / 15°
Telescopic distance of sheath	330mm
Patient Weight Capacity	150Kgs
Specification of guardrail	25*10mm /25*8mm /30*10mm /30*8mm /28.5*9.5mm
Product Weight	13kgs