



Tel : +86 15250358182
E-mail : info@medik-medical.com

ABS Bedside cabinet

YA-B07

Compact design, durable structure constructed in ABS material while with light weight, all the points make this bedside locker to be an ideal option for patient ward area at economical budget.



Castor

The bedside table has wheels for easy movement and two brakes. It can be fixed after moving to the designated position and is durable.

- Whole ABS forming
- Cabinet top, drawer, tray, layer board, door and the base are made of ABS
- Towel hangers on the 2 sides of cabinet, that can be adjustable
- Dinning board is between the drawer and tabletop
- Thermos , washbasin and other commodity can be put in it
- With four casters, easy to move

NO.	ITEM	PARAMETERS
01	Model	YA-B07
02	Size	L475*W475*H750mm
03	Accessories	1 drawer, 1 cabinet,1 sliding tray and shoe shelf
04	Caster	Dia.50mm


Material

The whole body is injection molded with 100% ABS plastic, which is non-toxic, strong, rust-proof, anti-corrosion, easy to clean, and has a long service life.



Sliding tray

With a hidden tray, there are grooves for placing water cups and chopsticks, which can increase the area of the table after sliding out.



Towel hook

Towel rail is with concealed design and can be stored when not in use, same two hanging hooks are with same design for usage and storage convenience.



ABS shoe holder

ABS shoe storage space which can put shoes in it.

