

# YA-D5-1

## Automatic Hospital ICU Bed



*Better care  
Better comfort!*

# YA-D5-1

# Automatic Hospital ICU Bed

## DESCRIPTION

MEDIK YA-D5-1 electric hospital bed is full electric operation type ICU bed to take care of the patient in ICU area. It will provide convenient support for both caregiver and patient. The road to recovery is often a critical and sensitive experience. Both patient and medical staff will have to work hand-in-hand to achieve their goal of nursing the patient back to health.

## FEATURES

- Bed Frame Electrophoresis and powder coating
- Bed platform removable ABS platform
- Motor Quiet and robust electric actuators provide reliable operation
- Build-in controller can control all electrically operated functions and movements
- 4 ABS side rails with Tuck-away model, collapsible below the bed surface
- Head and Footboard ABS material, detachable
- Central brake system with a wheel diameter of 125mm
- Hand remote control
- Patient control inner constructed in side rails
- Nurse control located in the footboard, with side rail control and a functions lockout feature
- Auto-counter to provide the patient from sliding onto foot end
- Electric CPR function control by one single button on the side rails
- Manual CPR

## STANDARD SPECIFICATION

Overall size	2200*980 *450-750mm
Lying size	1950*850 mm
Backrest adjustable	0-70°
Footrest adjustable	0-40°
Height adjustable	450-750mm
Trendelenburg	0-12°
Reverse-Trendelenburg	0-12°
Loading capacity	220kg



**BED BOARD**

The bed surface is detachable and has a handle for easy cleaning also it with vent holes for ventilation.



**SIDE RAIL LOCK**

Tuck-away full length railing meets IEC60601-2-52 safety standard, special protection mechanism can prevent patient to drop down the railing from inside bed.



**SIDERAIL CONTROL**

All functions can be operated by the built-in side rail control buttons, which are located inside and outside of the side rail, which is convenient for patients and medical staff to operate.



**Foot board control panel**

Nurse control panel embedded onto the foot board, railing control operation can be locked accordingly when position operation from patient is not wanted.



**Detachable head and foot board**

The detachable head and foot of the bed are equipped with safety locks to lock the head and foot board when it is not in use, and it can be quickly disassembled, which is convenient for medical staff to operate.



**CENTRAL LOCK**

Central lock system is convenient for the nurse to push the patient to the required position and quickly lock it with one key, which can appropriately reduce the workload of the medical staff.



## ABOUT US

Zhangjiagang Medi Medical Equipment Co.,Ltd, founded in 1995, with more than 20 years experience in developing hospital furniture, is one of the largest manufacturer of hospital bed, patient stretcher, medical cart and other related hospital furniture in China or even Asia.

Through scientific management complying with ISO9001&ISO13485, by adopting advanced facility like Germany TRUMPF CNC, Japan OTC Robot Welding, Auto-power coating line with Swiss GEMA gun , Medi focuses on providing brand quality product with CE approval.

Headquartered in Zhangjiagang nearby Shanghai, with subsidiary oversea offices in Hungary, Dubai, Thailand, New Zealand, Kenya, Colombia, Medi acting professionally and swiftly for both pre-sales and after-sales. The products are widely accepted by global clients, such as UK RIE hospital, UNIVERZITETNI KLINICNI CENTER, Kokilaben Dhirubhai Ambani Hospital, Singapore National University Hospital etc.



YOU TUBE



MEDIK

**Pls scan the QR code for more info**

**info@medik-medical.com | +86 199 6269 8550**

Medik reserves the right to make changes without notice in design, specifications and models. The only warranty Medik makes is the express written warranty extended on the sale or rental of its products.

© 2022 Medik Services, Inc. ALL RIGHTS RESERVED. APR53703 rev 1 01-JAN-2022 ENG – US