

# YA-D5-5

Electric Hospital Bed



*Better care  
Better comfort!*

# YA-D5-5

# Electric Hospital Bed

## DESCRIPTION

The YA-D5-5 bed has all the necessary functions of a standard ICU bed, ergonomic design, elegant appearance, and economical purchase in the hospital. The four-piece long guardrail design prevents the patient from falling and is safer.

## FEATURES

- Bed Frame Electrophoresis and powder coating
- Bed platform ABS material bed platform
- Motor Quiet and robust electric actuators provide reliable operation
- Hand remote controller of all electrically operated functions
- Handrails Four split side rail, drop down under the bed platform
- Detachable head/foot board with easy lock
- Castors Central brake system with wheel diameter of 125mm
- Auto-counter to provide the patient from sliding onto foot end
- Manual CPR on both sides(optional)
- Angle indicator With angle indicator embedded in the railing for back and Trendelenburg adjustment

## STANDARD SPECIFICATION

Overall size	2200*1060*510-800mm
Lying surface	1960*900mm
Height adjustable	510-800mm
Backrest adjustable	0° to 75°
Knee section adjustable	0° to 40°
Trendelenburg	0-12°
Reverse-trendeleburg	0-12°
Loading capacity	225kg



**CPR FUNCTION**

Manual CPR handles are accessible from both sides



**REMOTE CONTROL**

The hand control with clear icon is beneficial for medical staff which can control the bed position from a remote distance.



**BUMPER CASTOR**

There are bumper castors on the four corners of the bed, which can reduce the impact of accidentally hitting the wall during the movement.



**ANGLE INDICATOR**

Angle Indicator on the guardrail, backrest railing indicator for back section adjustment with visible 30 degree mark, knee railing indicator to show degree of trendelenburg



**CENTRAL LOCK**

Central lock system is convenient for the nurse to push the patient to the required position and quickly lock it with one key, which can appropriately reduce the workload of the medical staff.



**EMERGENCY "STOP" BUTTON**

The doctor can stop movements in case of emergency



**SIDE RAIL LOCK**

Tuck-away full length railing meets IEC60601-2-52 safety standard, special protection mechanism can prevent patient to drop down the railing from inside bed



**BED BOARD**

The ABS bed surface is detachable and has a handle for easy cleaning. The bed surface has vent holes for ventilation.

## ABOUT US

Zhangjiagang Medi Medical Equipment Co.,Ltd, founded in 1995, with more than 20 years experience in developing hospital furniture, is one of the largest manufacturer of hospital bed, patient stretcher, medical cart and other related hospital furniture in China or even Asia.

Through scientific management complying with ISO9001&ISO13485, by adopting advanced facility like Germany TRUMPF CNC, Japan OTC Robot Welding, Auto-power coating line with Swiss GEMA gun , Medi focuses on providing brand quality product with CE approval.

Headquartered in Zhangjiagang nearby Shanghai, with subsidiary oversea offices in Hungary, Dubai, Thailand, New Zealand, Kenya, Colombia, Medi acting professionally and swiftly for both pre-sales and after-sales. The products are widely accepted by global clients, such as UK RIE hospital, UNIVERZITETNI KLINICNI CENTER, Kokilaben Dhirubhai Ambani Hospital, Singapore National University Hospital etc.



YOU TUBE



MEDIK

**Pls scan the QR code for more info**

**info@medik-medical.com | +86 199 6269 8550**

Medik reserves the right to make changes without notice in design, specifications and models. The only warranty Medik makes is the express written warranty extended on the sale or rental of its products.

© 2022 Medik Services, Inc. ALL RIGHTS RESERVED. APR53703 rev 1 01-JAN-2022 ENG – US