

# YA-DH4-1

Rotating Bed



*Better care  
Better comfort!*

# YA-DH4-1

## Rotating Bed

### DESCRIPTION

YA-DH4-1 is a full electric control home care bed. It can provide comfortable position for aged person. The special rotation design will help aged person exit and return bed easily without care giver extra support. It is the reasonable solution for daily home using or home care center.

### FEATURES

- Cold rolled steel with powder coating treatment
- All function operation via hand remote control
- Detachable design for small and quick loading
- Wood material head/foot board
- Stable bearing design for chair position rotation

### STANDARD SPECIFICATION

Overall Size	2100*1000*510mm
Height adjustment range	510-810mm
Back board	80°
Knee board	45°
Chair position rotation	90 degree /One button control
Loading weight	240 kg
Packing weight	N.W 145 kg / G.W 165 kg



**Castors**

Dia 125mm swivel type castors, double side wheeling will provide heavy duty loading, 2 with locking.



**Hand remote controller**

All the function can be controlled via the remote hand control.



**Head/foot board**

Environment friendly MDF material, easy for daily clean.



**Side railing**

Side railing can be rotated 90 degree, it will provide more convenient support for aged person .

## ABOUT US

Zhangjiagang Medi Medical Equipment Co.,Ltd, founded in 1995, with more than 20 years experience in developing hospital furniture, is one of the largest manufacturer of hospital bed, patient stretcher, medical cart and other related hospital furniture in China or even Asia.

Through scientific management complying with ISO9001&ISO13485, by adopting advanced facility like Germany TRUMPF CNC, Japan OTC Robot Welding, Auto-power coating line with Swiss GEMA gun , Medi focuses on providing brand quality product with CE approval.

Headquartered in Zhangjiagang nearby Shanghai, with subsidiary oversea offices in Hungary, Dubai, Thailand, New Zealand, Kenya, Colombia, Medi acting professionally and swiftly for both pre-sales and after-sales. The products are widely accepted by global clients, such as UK RIE hospital, UNIVERZITETNI KLINICNI CENTER, Kokilaben Dhirubhai Ambani Hospital, Singapore National University Hospital etc.



YOU TUBE



MEDIK

**Pls scan the QR code for more info**

**info@medik-medical.com | +86 139 6223 1292**

Medik reserves the right to make changes without notice in design, specifications and models. The only warranty Medik makes is the express written warranty extended on the sale or rental of its products.

© 2022 Medik Services, Inc. ALL RIGHTS RESERVED. APR53703 rev 1 01-JAN-2022 ENG – US