

YA-DS-D04

Dialysis Chair



*Better care
Better comfort!*

YA-DS-D04

DIALYSIS CHAIR

DESCRIPTION

The YA-DS-D04 dialysis chair is ergonomically designed, which can better provide patients with more comfortable environment, reduce the psychological pressure of patients, and make the process of blood donation or dialysis more comfortable. Ideal for hemodialysis blood donation centers, chemotherapy and rehabilitation homes.

FEATURES

- Ergonomic design gives blood donors a sense of security and comfort and alleviate the psychological pressure, makes the whole process easier and enjoyable.
- The appearance is made of high-strength fiber material and poly-layer paint
- Upholstery with 60mm thick foam
- Chair frame is made of premium powder coated steel
- Four \varnothing 100 mm casters with individual locking system are included to facilitate the movement of the donor chair.

STANDARD SPECIFICATION

Overall size	1950*900*580-780mm
Seat size	520*600*580mm
Back section adjustment	-12°-75°
Leg section adjustment	-75°-12°
Height adjustment	580-780mm
Trendelenburg	0-12°
Height of armrest	230-290mm
Armrest swing out	0-180°
Pillow size	410*230*80mm
Loading capacity	240kg



Footrest

Foldable footrest can give foot support during a semi fowler position and prevent the patient from sliding down.



Cpr Handle

Manual CPR is behind the back section, which can quickly adjust the chair to anti-shock position for emergency treatment.



Hand Control

Equip with Denmark Linak brand motor control system, easily control via remote hand control.



Arm Rest

The cushioned armrest can move with the backboard when the back section is being adjusted, The upper armrest can swing out to facilitate the patient during treatment.



Pillow

The 8 cm thickness pillow provides head support during dialysis process, and it is also very convenient to remove for hygiene.



Castors

Four ø 100 mm casters with individual locking systems are included to facilitate the movement of the dialysis chair.



Iv Pole

The chair can be equipped with height-adjustable IV pole.



Dinner Table

The small dinner table is on the left side of the chair, which can be rotated 360°, specially designed with sundries and water cup.

ABOUT US

Zhangjiagang Medi Medical Equipment Co.,Ltd, founded in 1995, with more than 20 years experience in developing hospital furniture, is one of the largest manufacturer of hospital bed, patient stretcher, medical cart and other related hospital furniture in China or even Asia.

Through scientific management complying with ISO9001&ISO13485, by adopting advanced facility like Germany TRUMPF CNC, Japan OTC Robot Welding, Auto-power coating line with Swiss GEMA gun , Medi focuses on providing brand quality product with CE approval.

Headquartered in Zhangjiagang nearby Shanghai, with subsidiary oversea offices in Hungary, Dubai, Thailand, New Zealand, Kenya, Colombia, Medi acting professionally and swiftly for both pre-sales and after-sales.The products are widely accepted by global clients, such as UK RIE hospital, UNIVERZITETNI KLINICNI CENTER, Kokilaben Dhirubhai Ambani Hospital, Singapore National University Hospital etc.



YOU TUBE



MEDIK

Pls scan the QR code for more info

info@medik-medical.com | +86 139 6223 1292

Medik reserves the right to make changes without notice in design, specifications and models. The only warranty Medik makes is the express written warranty extended on the sale or rental of its products.

© 2022 Medik Services, Inc. ALL RIGHTS RESERVED. APR53703 rev 1 01-JAN-2022 ENG – US