

medik

Tel: +86 15250358182

E-mail: info@medik-medical.com

Zhangjiagang
Medi Medical Equipment
Co., Ltd

YA-DS-D06

/ Blood donor (dialysis) chair



YA-DS-D06 is a low-cost chair for blood donor or dialysis center. Electric operation can give patient 5 position choice during the treatment process. Ergonomic design provides patients with more comfortable experience. Ideal for dialysis centers, blood donation centers and nursing homes.



PUSH HAND

The push hand behind the back section can better enable the staff to move the chair.



PILLOW

High-density foam pillow can be easily attached and removable to the head section for comfort.



OVERBED TABLE (OPTIONAL)

The overbed table is optional, and necessary items can be placed on it and installed on the right side of the chair.



IV POLE HOLDER

High-quality epoxy-painted material is used to place the infusion stand, and there are knobs to fix it.



HAND CONTROL

Hand control with visible icons for easy operation. CPR button is incorporated for emergency case.



CASTOR

5" casters with individual brakes allow for better handling and control.

Technical parameter

Overall size	1000*880*1300mm
Back section adjustment	-10°~75°
Leg section adjustment	-80°~12°
Height adjustment	560-820mm
Trendelenburg position	0-12°
Loading capacity	200KG

Main features

- * Padded push-handle for easy maneuverability.
- * Extra-wide padded armrests maximize comfort
- * Durable epoxy coated welded steel frame is easy to clean and maintain
- * Equipped with 3 motors for height, back and leg rest adjustment
- * Multi position and inclination adjustments for the arm rests
- * Four 3 inch castors with individual locking system
- * Small dining table
- * With CPR function

