

# YA-DS-D07

Electric Dialysis Chair



*Better care  
Better comfort!*

# YA-DS-D07

# ELECTRIC DIALYSIS CHAIR

## DESCRIPTION

The YA-DS-D07 dialysis chair with ergonomic design to provide comfort during the dialysis treatment. The aptotic lying position will never bring abdominal pressure, it's an ideal luxury electric dialysis chair for hemodialysis centers, hospitals, retirement homes, etc.

## FEATURES

- Ergonomic design, no abdominal pressure
- Luxurious soft mattress to enhance patient comfort
- Trendelenburg, reverse Trendelenburg and any position in-between
- Adjust chair position via an electric knob
- Adjustable armrest
- Castors with brakes for easy movement

## TECHNICAL PARAMETER

|   |                |
|---|----------------|
| Model                                   | YA-DS-D07      |
| Overall size                            | 1900*580*550mm |
| Back section angle                      | 0°~65°         |
| Leg section angle                       | 0°~20°         |
| Armrest up and down adjustment distance | 0 ~ 100mm      |
| Armrest awing out angle                 | 0°~180°        |
| The maximum load weight                 | ≥ 180kg        |
| Castor                                  | Dia75mm        |

## PRODUCT DETAILS

### Trendelenburg position

Trendelenburg position adjusted by the knob, and can stop at any position.



### Recliner position

Recliner position adjusted by the knob.



### Reverse trendelenburg

Reverse trendelenburg position adjusted by the knob, and can stop at any position.



### Armrest

Height of armrest can be adjusted and can swing out.



## Knob adjustment

Chair position adjusted by the electric knob.



---

## Vinyl upholstery


High quality vinyl upholstery with antimicrobial/bacterial finish, allowing for durability and safety.



---

## Castors

4 castors of dia 75mm with individual brake, easy to move.



## ABOUT US

Zhangjiagang Medi Medical Equipment Co.,Ltd, founded in 1995, with more than 20 years experience in developing hospital furniture, is one of the largest manufacturer of hospital bed, patient stretcher, medical cart and other related hospital furniture in China or even Asia.

Through scientific management complying with ISO9001&ISO13485, by adopting advanced facility like Germany TRUMPF CNC, Japan OTC Robot Welding, Auto-power coating line with Swiss GEMA gun , Medi focuses on providing brand quality product with CE approval.

Headquartered in Zhangjiagang nearby Shanghai, with subsidiary oversea offices in Hungary, Dubai, Thailand, New Zealand, Kenya, Colombia, Medi acting professionally and swiftly for both pre-sales and after-sales.The products are widely accepted by global clients, such as UK RIE hospital, UNIVERZITETNI KLINICNI CENTER, Kokilaben Dhirubhai Ambani Hospital, Singapore National University Hospital etc.



YOU TUBE



MEDIK

**Pls scan the QR code for more info**

**info@medik-medical.com | +86 139 6223 1292**

Medik reserves the right to make changes without notice in design, specifications and models. The only warranty Medik makes is the express written warranty extended on the sale or rental of its products.