

# PRODUCT FUNCTION

YA-DS-D09 Blood donor (dialysis) provides dialysis center a comfortable product at very competitive cost. Ergonomic and compact design not only provides comfort for patient but also save the space for dialysis center. Motorized operation can give patient more position choice during the treatment process. Ideal for long term care facilities, retirement homes and nursing homes.



Fowler position



Trendelenburg position



Full Flat Position



## YA-DS-D09

DIALYSIS CHAIR

Zhangjiagang  
Medi Medical Equipment Co.,Ltd

[www.medik-medical.com](http://www.medik-medical.com)



You Tube Video



Medik Company

Address : No.3 Jinxiu Road,Jinfeng Industrial Zone,Jinfeng Town,Zhangjiagang,Jiangsu Province,China  
Tel : +86 15250358182 / Fax : +86 512 58688450  
E-mail : [info@medik-medical.com](mailto:info@medik-medical.com)

BETTER CARE | BETTER COMFORT

# PRODUCT TECHNICAL PARAMETERS

Medi Medical Equipment Co., Ltd  
Better care Better comfort

## PARAMETER LIST

Here are the main parameters of the product function

NO.	ITEM	PARAMETER
1	Model	YA-DS-D09
2	Name	Dialysis chair
3	Overall size	920*850*560mm
4	Seat size	570*500mm
5	Back section adjustment	-10°~75°
6	Leg section adjustment	-80°~0°
7	Trendelenburg position	0-12°
8	Loading capacity	200 KG

BETTER CARE BETTER COMFORT



## PRODUCT DETAILS



### Pillow

Removable pillow can be easily attached to head section for comfort.



### Armrest

The cushioned armrest can move with the backboard when the back-section is being adjusted



### Hand control

Hand control with visible icons for easy operation. CPR button is incorporated for emergency case.



### Footrest

Foldable footrest can give foot support during a semi fowler position and prevent the patient from sliding downwards.



### Casters

Elegant design silent durable caster with individual brake



### IV Pole holder

High-quality epoxy-painted material is used to place the infusion stand, and there are knobs to fix it. IV pole holder diameter is 25mm