

YA-DS-M04D

Phlebotomy Chair



*Better care
Better comfort!*

YA-DS-M04D PHLEBOTOMY CHAIR

DESCRIPTION

YA-DS-M04D is a functionally designed medical chair with a storage drawer and adjustable arm rest which is widely used in different hospital fields such as blood donation and vaccination. The ergonomic design provides a comfortable experience for patients and nurses.

FEATURES

- Powder coated steel structure
 - Back and Seat made of foam upholstered with PVC
 - Padded armrest adjustable 90 degree
 - With a drawer at side for storage
 - Knock down construction
-

STANDARD SPECIFICATION

Overall size	680*650*950mm
Seat size	550*520*560mm
Loading capacity	180kg



Arm rest adjust

The arm rest can be slid according to the needs of the users and flip up to facilitate leaving the seat.



Arm support adjust

The arm section can be adjust by the knob, with different height to provide more comfortable experience to the users.



Storage drawer

The epoxy powder coated drawers can store items, and the drawer is padded to comfortably place your hands on it.



Leg protection

Four legs with rubber feet to protect and prevent slipping.



Material

Back and Seat are made of foam upholstered with PVC.

ABOUT US

Zhangjiagang Medi Medical Equipment Co.,Ltd, founded in 1995, with more than 20 years experience in developing hospital furniture, is one of the largest manufacturer of hospital bed, patient stretcher, medical cart and other related hospital furniture in China or even Asia.

Through scientific management complying with ISO9001&ISO13485, by adopting advanced facility like Germany TRUMPF CNC, Japan OTC Robot Welding, Auto-power coating line with Swiss GEMA gun , Medi focuses on providing brand quality product with CE approval.

Headquartered in Zhangjiagang nearby Shanghai, with subsidiary oversea offices in Hungary, Dubai, Thailand, New Zealand, Kenya, Colombia, Medi acting professionally and swiftly for both pre-sales and after-sales.The products are widely accepted by global clients, such as UK RIE hospital, UNIVERZITETNI KLINICNI CENTER, Kokilaben Dhirubhai Ambani Hospital, Singapore National University Hospital etc.



YOU TUBE



MEDIK

Pls scan the QR code for more info

info@medik-medical.com | +86 139 6223 1292

Medik reserves the right to make changes without notice in design, specifications and models. The only warranty Medik makes is the express written warranty extended on the sale or rental of its products.

© 2022 Medik Services, Inc. ALL RIGHTS RESERVED. APR53703 rev 1 01-JAN-2022 ENG – US