






# YA-DS-R02

RECLINING PHLEBOTOMY CHAIR



## PRODUCT DETAILS



### Head rest

Comfortable built-in headrest



### Casters

5" Total lock casters allow for better handling and control.



### Upholstered arm supports

Upholstered arm supports provide extra comfort in recline position



### Non-porous trays

Non-porous trays, stores on side of chair.



### Comfortable armrests

Comfortable armrests are angled to provide support in any reclined position



### Release handle

Gas cylinder release handle for quick back release.

- Four positions: upright, TV, full recline and Trendelenburg
- Enhanced contouring and cushioning ideal for long-term seating and care
- Gas cylinder controls deep recline and Trendelenburg
- Side panels can be removed for quick, thorough cleaning
- Non-porous trays, stores on side of chair.
- Gas cylinder release handle for quick back release
- Comfortable armrests are angled to provide support in any reclined position
- Sturdy, easy-to-clean steel frame, 113kg weight capacity
- 5" Casters (2 with locks)

1	Model	YA-DS-R02
2	Overall size (H*W*D)	1205 x 775 x 965 mm
3	Seat size (H*W*D)	545 x 545 x 508 mm
4	Floor to Seat Height	675mm
5	Width between arms	560mm
6	Loading capacity	113kg
7	Product weight	53kg

