

medik

Tel: +86 15250358182

E-mail: info@medik-medical.com

Zhangjiagang
Medi Medical Equipment
Co., Ltd

YA-EC-M03

/ Exam couch

Exam couch YA-EC-M03 has sturdy sprayed frame, equipped with dirt basin and drawers for convenient storage of items for medical staff and is suitable for general inspection, endoscopy, etc.



BACK SECTION

The back section is adjusted by ratchet system to be between 0~70°



DRAWERS

Exam couch is equipped with 2 drawers with central lock



ROLL PAPER HOLDER

Roll paper holder is installed at the side of head section for easy installation and replacement of roll paper



LATHER

leather is made of medical artificial leather with flame-retardant, water-proof and wear-resistant



BASIN

Dirt basin under the back plate, which can be pulled out

Main features

- * Manufactured in epoxy-painted steel
- * Top is upholstered and covered with washable plastic material
- * Complete with 2 drawers, paper roll holder and s.s. tray stainless steel basin for fluid collection under the back
- * Mounted on Ø 80 mm.wheels with brake

Technical parameter

Overall size	1800*800*550 mm
Back section adjustment	0~75°
Loading capacity	150kg

