






**MAIN FEATURES**

- \* Electrophoresis and powder coating
- \* Head & Foot board:ABS,detachable
- \* Side rails:2pcs,Aluminum alloy,collapsible
- \* Castor:diameter of 125mm,two with brakes
- \* I.V.pole



The YA-M2-3 two cranks manual hospital bed is a manual crank system offering function of back and knee up & down. More over the medical adjustable bed is independent castor locking system.

**PRODUCT DETAILS**



**TECHNICAL PARAMETER**

ITEM	PARAMETER
Model	YA-M2-3
Overall size	2100*980*600mm
Lying size	1910*830mm
Back section up/down	0-75°
Foot section up/down	0-40°
Fixed height	600mm
Safe loading capacity	200kg