

medik

Tel : +86 15250358182
E-mail : info@medik-medical.com




YA-M3-1

MANUAL HOSPITAL BED




MAIN FEATURES

- * Bed Frame electrophoresis and powder coating
- * Bed platform MS perforated platform with ventilation hole, molded
- * Handrails collapsible aluminum railing with 6 bars
- * Head and Foot board detachable head/foot board with easy lock
- * Cranks 3 peaces foldable crank
- * Castors Caster diameter of 125mm with diagoal locking
- * Bumper castors bumper castors on 4 corners of bed



BACK SECTION



FOOT SECTION



HEIGHT ADJUSTMENT

PRODUCT DETAILS



FOLDABLE CRANK HANDLE



CASTORS



HEADBOARD / FOOTBOARD



BUMPER COVER



IV POLE



STEEL | PAINTING

TECHNICAL PARAMETER

ITEM	PARAMETER
Model	YA-M3-1
Overall size	2150*980*450-750mm
Back section adjustable	0-75°
Foot section adjustable	0-40°
High-low	450-750mm
Lying size	1900*860mm
Safe working load	200kg