

# YA-PD5-5

Pediatric Bed



*Better care  
Better comfort!*

YA-PD5-5

Pediatric Bed

DESCRIPTION

The YA-PD5-5 is a pediatric Bed with electric and manual controls. The extra height fence design provides the children with great safety, up-down sliding rails with double lock prevents children from opening the rail accidentally which is another key safety design.

FEATURES

- Hand remote control
- Head and footboard: Powder spray
- Side rails: Two powder spray side rails with HI-LOW position
- Castors: Four 125mm castors and two with brakes
- Mattress is also available
- Backrest & footrest adjust by the manual crank
- With I.V. pole

STANDARD  
SPECIFICATION

Model	YA-PD5-5
Overall size	1300*650*1550mm
Height adjustable	450-800mm
Backrest adjustable	0° to 75°
Footrest adjustable	0° to 40°
Trendelenburg	0° to 12°
Reverse trendelenburg	0° to 12°
Loading capacity	120kg



### **SIDE RAIL**

The side rail with 3 gears to adjust the height, anti-pinch hands to protect children's safety . When adjusting, It needs to be pulled two handles at the same time to adjust the side rail which prevents children from misoperation.



### **REMOTE CONTROL**

Height and Trendelenburg & Reverse-Trendelenburg functions are electrically controlled. The icons on the remote control are clear and easy for nurses to operate.



### **CASTORS**

125mm diameter castors easy for nurses to push the pediatric bed to any position they want , two with brakes.



### **MANUAL CRANK**

Backrest & footrest functions are adjusted by the manual crank .



## ABOUT US

Zhangjiagang Medi Medical Equipment Co.,Ltd, founded in 1995, with more than 20 years experience in developing hospital furniture, is one of the largest manufacturer of hospital bed, patient stretcher, medical cart and other related hospital furniture in China or even Asia.

Through scientific management complying with ISO9001&ISO13485, by adopting advanced facility like Germany TRUMPF CNC, Japan OTC Robot Welding, Auto-power coating line with Swiss GEMA gun , Medi focuses on providing brand quality product with CE approval.

Headquartered in Zhangjiagang nearby Shanghai, with subsidiary oversea offices in Hungary, Dubai, Thailand, New Zealand, Kenya, Colombia, Medi acting professionally and swiftly for both pre-sales and after-sales. The products are widely accepted by global clients, such as UK RIE hospital, UNIVERZITETNI KLINICNI CENTER, Kokilaben Dhirubhai Ambani Hospital, Singapore National University Hospital etc.



YOU TUBE



MEDIK

**Pls scan the QR code for more info**

**info@medik-medical.com | +86 199 6269 8550**

Medik reserves the right to make changes without notice in design, specifications and models. The only warranty Medik makes is the express written warranty extended on the sale or rental of its products.

© 2022 Medik Services, Inc. ALL RIGHTS RESERVED. APR53703 rev 1 01-JAN-2022 ENG – US