

YA-PM3-2

Pediatric Bed



*Better care
Better comfort!*

YA-PM3-2

Pediatric Bed

DESCRIPTION

The YA-PM3-2 bears all the necessary features for a standard Pediatric Bed, and the high guardrail is designed to prevent children from falling from the bedside. Comfortably for children' recovery, economically for hospital to buy.

FEATURES

- Powder Spray Head and footboard
- Side rails: Two powder spray side rails with HI-LOW position
- Castors: Four 100mm castors and two with brakes
- With I.V. pole
- Mattress is also available

STANDARD SPECIFICATION

Model	YA-PM3-2
Overall size	1300*650*1550mm
Height adjustable	450-800mm
Backrest adjustable	0° to 75°
Trendelenburg	0° to 12°
Reverse trendelenburg	0° to 12°
Loading capacity	120kg



SIDE RAIL PROTECTION MECHANISM

The overall height of the guardrail is adjustable, and the raised guardrail protection mechanism prevents children from pushing the guardrail and preventing them from falling off the bed.



SIDE RAIL

The height of the Side rail can prevent children from standing up and falling, the gap is small to avoid falling.



CASTORS

100 mm diameter castors easy for nurses to push the pediatric bed to any position they want , two with brakes.



BED BOARD

The welded bed board structure is more breathable and convenient for nurses to clean.



ABOUT US

Zhangjiagang Medi Medical Equipment Co.,ltd, founded in 1995, with more than 20 years experience in developing hospital furniture, is one of the largest manufacturer of hospital bed, patient stretcher, medical cart and other related hospital furniture in China or even Asia.

Through scientific management complying with ISO9001&ISO13485, by adopting advanced facility like Germany TRUMPF CNC, Japan OTC Robot Welding, Auto-power coating line with Swiss GEMA gun , Medi focuses on providing brand quality product with CE approval.

Headquartered in Zhangjiagang nearby Shanghai, with subsidiary oversea offices in Hungary, Dubai, Thailand, New Zealand, Kenya, Colombia, Medi acting professionally and swiftly for both pre-sales and after-sales.The products are widely accepted by global clients, such as UK RIE hospital, UNIVERZITETNI KLINICNI CENTER, Kokilaben Dhirubhai Ambani Hospital, Singapore National University Hospital etc.



YOU TUBE



MEDIK

Pls scan the QR code for more info

info@medik-medical.com | +86 199 6269 8550

Medik reserves the right to make changes without notice in design, specifications and models. The only warranty Medik makes is the express written warranty extended on the sale or rental of its products.

© 2022 Medik Services, Inc. ALL RIGHTS RESERVED. APR53703 rev 1 01-JAN-2022 ENG – US