



### Back section

The back section is adjusted by gas spring to be between 0-75°



### Foot pedal

Height adjustment and Trendelenburg and reverse-Trendelenburg function are controlled by the foot pedal, which can release the hands of the medical staff



### Handle

Ergonomic handle located at the head end can let the staff to transport the stretcher in any way.



### Oxygen holder

oxygen cylinder holder is specially designed on the base of the stretcher



### Caster

Diameter 200mm castor with central brake system is excellent for transportation



### 5th direction wheel

The stretcher has 5th direction wheel provides the stretcher with ideal manipulation and stability



### Restraint belts

The stretcher is equipped with two straps which are fixed on the bed surface by screws, which can well fix the patient and prevent the patient from shaking during the movement.



### IV pole hole

IV pole holes are installed on both sides of the stretcher

# YA-PS04 / PATIENT STRETCHER

YA-PS04 is a standard hydraulic patient stretcher designed for patient comfort and caregiver convenience which meet the needs of any healthcare facility. It's ergonomically positioned handles make it easy to maneuver.

**medik**

Tel : +86 15250358182  
E-mail : info@medik-medical.com



## Technical parameter

Overall size (L*W*H)	193*70*63-93cm
Back section adjustment	0-75°
Caster (Diameter)	200mm
Height adjustment	63-93cm
Trendelenburg	0-15°
Reverse Trendelenburg	0-15°
Loading capacity	180 kg

## Function

- \* Back section
- \* HI-LOW
- \* Reverse Trendelenburg
- \* Trendelenburg

## Main features

- \* Metal parts coated with electrostatic epoxy polyester powder
- \* IV pole and protective bumpers on the all corners
- \* 200mm diameter castors capable of turning 360°
- \* Trendelenburg and reverse Trendelenburg are controlled with hydraulic pump
- \* Designed ABS plastic base cover is dust proof, and can be placed oxygen cylinders
- \* Body fixing strap can prevent the patient from shaking during the movement

