



## MAIN FEATURES

YA-S02 Doctor stool is designed to provide a convenient and dependable solution at affordable cost. The Stool features a rugged, five-leg base with swivel casters that provide a generous weight capacity of 300 lbs.

- \* PU with metal insert
- \* Five-Leg Base with Swivel Casters for Excellent Mobility
- \* High Density Foam Padding
- \* Rugged Vinyl Top Upholstery

**medik**  
Tel : +86 15250358182  
E-mail : info@medik-medical.com



**YA-S02**  
DOCTOR STOOL



### Height Adjustment

Height adjustable is from 430mm to 530mm

①



### Base

Five-leg base strong and stable

②



### Castor

Dia 50mm castors move smoothly

③



### Seat Surface

PU leather, soft and comfortable

④



## TECHNICAL PARAMETER

ITEM	PARAMETER
Seat	Dia 350mm*H10mm
Gas lift	100mm travel distance
Base	Dia 550mm
Castor	Dia 50mm