

YA-S08

Doctor Stool



*Better care
Better comfort!*

YA-S08

DOCTOR STOOL

DESCRIPTION

With a stable disc-shaped base, easy gas spring height adjustment, and thick seat cushions and backrests, the YA-S08 nurse stool combines ergonomic comfort design with a beautiful appearance, making it suitable for a variety of scenarios in clinic and hospital.

FEATURES

- Thick foam padding with PU cover, antibacterial and easy to clean.
- Ergonomic design, with backrest.
- Height adjusted by gas spring.
- Stable large diameter disc base with swivel castors.

STANDARD SPECIFICATION

Overall size (from ground to seat surface)	Ø55*55-65cm
Height range (from ground to seat surface)	55~65cm
Back rest height from seat surface	28cm
Seat top diameter	Dia 400mm
Bottom base diameter	Dia 550mm
Castor	Dia 50mm



Height Adjustment

Height adjusted by gas spring,easy and durable.



Seat Cushions

Thick foam padding with PU cover ,antibacterial and easy to clean.



Backrest

Backrest with thick foam padding,removable.



Base And Castors

4 swivel castors of dia 50mm with big round base ,durable and easy to move.

ABOUT US

Zhangjiagang Medi Medical Equipment Co.,Ltd, founded in 1995, with more than 20 years experience in developing hospital furniture, is one of the largest manufacturer of hospital bed, patient stretcher, medical cart and other related hospital furniture in China or even Asia.

Through scientific management complying with ISO9001&ISO13485, by adopting advanced facility like Germany TRUMPF CNC, Japan OTC Robot Welding, Auto-power coating line with Swiss GEMA gun , Medi focuses on providing brand quality product with CE approval.

Headquartered in Zhangjiagang nearby Shanghai, with subsidiary oversea offices in Hungary, Dubai, Thailand, New Zealand, Kenya, Colombia, Medi acting professionally and swiftly for both pre-sales and after-sales.The products are widely accepted by global clients, such as UK RIE hospital, UNIVERZITETNI KLINICNI CENTER, Kokilaben Dhirubhai Ambani Hospital, Singapore National University Hospital etc.



YOU TUBE



MEDIK

Pls scan the QR code for more info

info@medik-medical.com | +86 139 6223 1292

Medik reserves the right to make changes without notice in design, specifications and models. The only warranty Medik makes is the express written warranty extended on the sale or rental of its products.

© 2022 Medik Services, Inc. ALL RIGHTS RESERVED. APR53703 rev 1 01-JAN-2022 ENG – US