

YA-SS08

Stair stretcher



*Better care
Better comfort!*

YA-SS08

STAIR STRETCHER

DESCRIPTION

To transfer the patient or wounded person is a big challenge when there is no lift in the building. While the Stair Stretcher YA-SS08 with motorized belt is a creative assistant to help the single operator to achieve the transferring task effortlessly and efficiently on staircases.

FEATURES

- With a compact and foldable structure, ensuring ease of storage and portability
- Secure and adjustable straps to offer a snug fit without compromising on comfort
- Incorporates a non-slip track to enhance patient safety during transport
- Compact and foldable structure, ensuring ease of storage and portability.
- Crafted from durable materials, the stretcher guarantees robustness and reliability during transportation.
- The Engineered frame to withstand the rigors of stair-related transport

TECHNICAL PARAMETER

Model	YA-SS08
Size (Unfolded)	1250*540*1580mm
Size (Folded)	1050*540*280mm
Loading capacity	159kg
Carton Size (1pc/carton)	1200*580*350mm



Adjustable tracks

Effortlessly tackle staircases with our stretcher's adjustable track angles for smooth and adaptable transfers.



Castors

The four wheels offer superior mobility, ensuring efficient and secure transfers in diverse conditions



Retractable handle

Facilitate smooth and convenient transfers with retractable handle, designed for easy handling.



Foldable arm rest

Improve caregiver efficiency with a foldable armrest, delivering ergonomic support in diverse scenarios, enhance user experience.



Foldable design

Prioritize space economy with the overall foldable design, allowing for efficient storage and quick deployment, ideal for emergency state of affairs.

ABOUT US

Zhangjiagang Medi Medical Equipment Co.,Ltd, founded in 1995, with more than 20 years experience in developing hospital furniture, is one of the largest manufacturer of hospital bed, patient stretcher, medical cart and other related hospital furniture in China or even Asia.

Through scientific management complying with ISO9001&ISO13485, by adopting advanced facility like Germany TRUMPF CNC, Japan OTC Robot Welding, Auto-power coating line with Swiss GEMA gun , Medi focuses on providing brand quality product with CE approval.

Headquartered in Zhangjiagang nearby Shanghai, with subsidiary oversea offices in Hungary, Dubai, Thailand, New Zealand, Kenya, Colombia, Medi acting professionally and swiftly for both pre-sales and after-sales.The products are widely accepted by global clients, such as UK RIE hospital, UNIVERZITETNI KLINICNI CENTER, Kokilaben Dhirubhai Ambani Hospital, Singapore National University Hospital etc.



YOU TUBE



MEDIK

Pls scan the QR code for more info

info@medik-medical.com | +86 139 6223 1292

Medik reserves the right to make changes without notice in design, specifications and models. The only warranty Medik makes is the express written warranty extended on the sale or rental of its products.